

Tube	A	B	C	D	E	F	G	H	I	J	K	L	M
Desired pH	4.10	4.20	4.30	4.40	4.50	4.60	4.70	4.80	4.90	5.00	5.10	5.20	5.30
Actual pH	_____	_____	_____	_____	_____	_____	_____	_____	_____	_____	_____	_____	_____
Added invertase at (<i>hh:mm:ss</i>) Vortex briefly; Move to 55 °C water bath.	_____	_____	_____	_____	_____	_____	_____	_____	_____	_____	_____	_____	_____
Added 3,5-DNS at (<i>hh:mm:ss</i>) Move to 100 °C water bath. Wait 5 minutes or so. Let tube cool.	_____	_____	_____	_____	_____	_____	_____	_____	_____	_____	_____	_____	_____
Absorbance at 540 nm													
Tube 1	_____	_____	_____	_____	_____	_____	_____	_____	_____	_____	_____	_____	_____
Tube 2	_____	_____	_____	_____	_____	_____	_____	_____	_____	_____	_____	_____	_____